



Holloware

Healthcare Excellence, Every Step of the Way.

Tahir Group Of Industries



+92-52-3257151



+92-300-666-9737



www.tgind.com



TAHIR GROUP OF INDUSTRIES



TAHIR GROUP OF INDUSTRIES



TAHIR GROUP OF INDUSTRIES



info@tgind.com
tgindsale@gmail.com



New Miana Pura East, Roras
Road P.O.S.I.E Sialkot
51310 Pakistan



CERTIFICATE

No. 920229



This is to certify that the Quality Management System for Medical Devices of the company



www.tgind.com

Tahir Group of Industries

Dubai Islamic Bank Street Opposite S.I.E, Shahab Pura Road, Shahabpura, Sialkot,
Punjab 51310
Pakistan

has been assessed and found to be in compliance with the Standard

EN ISO 13485:2016

applicable to

**Manufacturer and Exporter of Surgical, Dental,
Electrosurgical, Orthopedic, Laparoscopic and
Diagnostic Instruments, Halloware**

The certificate has been issued under No. **920229** for the registration period from
05 December 2024 to 04 December 2027.
The first certificate date of issue is 29 March 2023.


Approved by


Printed by



validity code: **D60F2786-742**

Check the validity of this certificate using this code at www.ll-c.info

CERTIFICATE

No. 920229

This is to certify that the Quality Management System of the company



www.tgind.com

Tahir Group of Industries

Dubai Islamic Bank Street Opposite S.I.E, Shahab Pura Road, Shahabpura, Sialkot,
Punjab 51310
Pakistan

has been assessed and found to be in compliance with the Standard

EN ISO 9001:2015

applicable to

**Manufacturer and Exporter of Surgical, Dental,
Electrosurgical, Orthopedic, Laparoscopic, Diagnostic
and Beauty Instruments, Halloware**

The certificate has been issued under No. **920229** for the registration period from
05 December 2024 to 04 December 2027.
The first certificate date of issue is 29 March 2023.


Approved by


Printed by



validity code: **A1130737-011**

Check the validity of this certificate using this code at www.ll-c.info

CERTIFICATE OF COMPLIANCE

This is to certify, that the hereunder described items of Medical Devices Directive 93/42/EEC as amended by Directive 2007/42/EC have proven their conformity to the Safety and Health requirements.

Manufactured by: **TAHIR GROUP OF INDUSTRIES**

Address: Dubai Islamic Bank Street Opposite S.I.E, Shahab Pura Road, Shahabpura, Sialkot-51310, Pakistan.

Product: **Surgical, Dental, Electrosurgical, Orthopedic, Laparoscopic, Diagnostic, Veterinary and Beauty Instruments, Mackintosh Rubber Sheet & Medical Silicone & Rubber Products.**

Classification: Class I
non-sterile, non-powered and non-measuring devices
(according, to the Manufacturer's declaration)

The product liability rests with the manufacturer under applicable regulations and standards, after fulfilling the relevant EU legislation requirements and withdrawing an EC declaration of conformity as per Annexure-VII, module A of the product, you are therefore licensed to CE mark the product(s) listed above per Article 17 of the Medical Device Directive. The manufacturer is responsible for the CE Marking process and is not exempted from carrying out all necessary compliance activities.

This certificate shall not be reproduced except in full and remains the property of CNC Services which must be returned on request. The holder must inform CNC Services, of any substantial changes that occurred in the product or process in order to examine whether this certificate remains valid.



For CNC Services



[Signature]
Project Engineer

Date of Issue: 24th March 2024
Date of Expiry: 23rd March 2026
Certificate Number: EC/5231/23

Welcome to Tahir Group Of Industries Hospital Hollowares



is a Hospital & Medical Holloware Supplier in The Worldwide, Committed to Providing Customers With The Highest Quality Products. Our Top Priority is Patient Safety, so we Constantly Ensure Customer Satisfaction and Deliver The Most Modern Technology . Tahir Group Of Industries Use High Grade Stainless Steel 304 & 1050 Quality H24 Aluminium.130° to 135° Temperature Sterilize

INDEX

- Kidney Trays and Wash Basin
- Bowls and Medicine cup
- Instruments Trays
- Instruments Trays, Boxes
- PRF Box, Endo Box Dental Cassette
- Mesh Trays and Cassette
- Surgical Box, Microsurgery Box, Mats
- Needle Cases, Petri Dishes, Cotton Jars
- Cotton Dressing Drums
- Jugs, Feeding & Septum Cup, Strainer
- Bed Pans and Urinal Pot
- Hospital Bucket & Urinal, Can



Kidney Trays and Wash Basin

Kidney / Emesis Trays

English pattern

Without cover	Length x Width x Height
TG-HW-0001	6" x 3½" x 1¾" (150 x 90 x 45 mm)
TG-HW-0002	8" x 3½" x 1¾" (200 x 95 x 45 mm)
TG-HW-0003	10" x 4½" x 2" (255 x 115 x 50 mm)



Kidney / Emesis Trays

German pattern

Without cover	Length x Height
TG-HW-0004	6" x 1¼" (150 x 30 mm)
TG-HW-0005	6¾" x 1½" (170 x 35 mm)
TG-HW-0006	8" x 1½" (200 x 35 mm)
TG-HW-0007	10" x 1¾" (250 x 40 mm)
TG-HW-0008	10¾" x 1½" (275 x 45 mm)
TG-HW-0009	12" x 2" (300 x 50 mm)
TG-HW-00010	16" x 2¾" (400 x 60 mm)



Gallipots

Without cover	With cover	Diameter x Height	OZ
TG-HW-00011	TG-HW-00023	1¾" x 1¼" (42 x 30 mm)	1
TG-HW-00012	TG-HW-00024	2" x 1½" (53 x 34 mm)	2
TG-HW-00013	TG-HW-00025	2¼" x 1¾" (58 x 40 mm)	3
TG-HW-00014	TG-HW-00026	2½" x 1¾" (64 x 48 mm)	4
TG-HW-00015	TG-HW-00027	2½" x 1¾" (68 x 47 mm)	5
TG-HW-00016	TG-HW-00028	2¾" x 1¾" (71 x 46 mm)	6
TG-HW-00017	TG-HW-00029	2¾" x 2" (71 x 50 mm)	7
TG-HW-00018	TG-HW-00030	3¼" x 2" (81 x 51 mm)	8
TG-HW-00019	TG-HW-00031	3¼" x 1½" (81 x 54 mm)	9
TG-HW-00020	TG-HW-00032	3½" x 1¾" (89 x 47 mm)	10
TG-HW-00021	TG-HW-00033	4" x 1¾" (100 x 48 mm)	11
TG-HW-00022	TG-HW-00034	4" x 1½" (100 x 53 mm)	12



Iodine Cups

	Diameter x Height	OZ
TG-HW-00035	3½" x 2" (90 x 50 mm)	2
TG-HW-00036	(90 x 50 mm)	6
TG-HW-00037	(90 x 50 mm)	14



Medicine Cups

Graduated

	Diameter x Height	ml
TG-HW-00037	2¼" x 1¼" (58 x 30 mm)	25
TG-HW-00038	2" x 2" (52 x 50 mm)	50



Bowls and Medicine cup

Lotion Bowls

Without cover	With cover	Diameter x Height
TG-HW-00039	TG-HW-00047	1½" x 3" (45 x 75 mm)
TG-HW-00040	TG-HW-00048	1½" x 4" (50 x 100 mm)
TG-HW-00041	TG-HW-00049	2½" x 5" (60 x 125 mm)
TG-HW-00042	TG-HW-00050	2¾" x 6" (65 x 150 mm)
TG-HW-00043	TG-HW-00051	2¾" x 7" (70 x 180 mm)
TG-HW-00044	TG-HW-00052	3" x 8" (75 x 200 mm)
TG-HW-00045	TG-HW-00053	4" x 6" (100 x 150 mm)
TG-HW-00046	TG-HW-00054	4½" x 12" (110 x 300 mm)



Sponge Bowls

	Diameter x Height
TG-HW-00055	4½" x 2" (115 x 50 mm)
TG-HW-00056	4¾" x 2¼" (120 x 55 mm)
TG-HW-00057	5¾" x 2" (145 x 50 mm)
TG-HW-00058	6" x 2½" (150 x 65 mm)



Solution Bowls

	Diameter x Height
TG-HW-00059	7" x 3" (18 x 75 mm)
TG-HW-00060	6½" x 3" (165 x 75 mm)
TG-HW-00061	7½" x 3" (190 x 75 mm)
TG-HW-00062	8" x 3¼" (200 x 8 mm)
TG-HW-00063	9¾" x 4" (245 x 100 mm)
TG-HW-00064	10" x 3½" (250 x 90 mm)
TG-HW-00065	12" x 4½" (300 x 110 mm)
TG-HW-00066	12¾" x 5¼" (325 x 135 mm)



Solution Bowls

Tapered

	Diameter x Height
TG-HW-00067	6½" x 2½" (160 x 60 mm)
TG-HW-00068	7" x 2¾" (175 x 70 mm)
TG-HW-00069	8½" x 3¼" (220 x 80 mm)
TG-HW-00070	8½" x 3½" (225 x 90 mm)
TG-HW-00071	9½" x 4½" (240 x 110 mm)
TG-HW-00072	10" x 4½" (255 x 115 mm)
TG-HW-00073	10¾" x 4½" (275 x 115 mm)
TG-HW-00074	11½" x 5" (290 x 125 mm)
TG-HW-00075	12½" x 4¾" (310 x 120 mm)
TG-HW-00076	13¾" x 5" (345 x 125 mm)
TG-HW-00077	15¼" x 5" (380 x 130 mm)
TG-HW-00078	16½" x 5½" (415 x 140 mm)



Wash Basins

	Diameter x Height
TG-HW-00079	12" x 4" (300 x 100 mm)
TG-HW-00080	14" x 4" (350 x 100 mm)
TG-HW-00081	16" x 4" (400 x 100 mm)



Instruments Trays

Instruments Trays

Without cover	With cover	Diameter x Height	
TG-HW-00082	TG-HW-000100	6" x 4" x 2"	(150 x 100 x 50 mm)
TG-HW-00083	TG-HW-000101	8" x 6" x 2"	(200 x 150 x 50 mm)
TG-HW-00084	TG-HW-000102	9" x 6" x 1 1/4"	(230 x 155 x 30 mm)
TG-HW-00085	TG-HW-000103	10" x 8" x 2"	(250 x 200 x 50 mm)
TG-HW-00086	TG-HW-000104	12" x 6" x 2"	(300 x 150 x 50 mm)
TG-HW-00087	TG-HW-000105	12" x 7" x 1 1/4"	(300 x 175 x 30 mm)
TG-HW-00088	TG-HW-000106	12" x 7" x 1 3/4"	(300 x 175 x 40 mm)
TG-HW-00089	TG-HW-000107	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-00090	TG-HW-000108	12" x 10" x 2"	(300 x 250 x 50 mm)
TG-HW-00091	TG-HW-000109	14" x 8 1/2" x 2"	(350 x 220 x 50 mm)
TG-HW-00092	TG-HW-000110	14" x 10" x 1 1/2"	(350 x 250 x 36 mm)
TG-HW-00093	TG-HW-000111	14 1/2" x 12" x 2"	(355 x 300 x 50 mm)
TG-HW-00094	TG-HW-000112	16" x 12" x 2"	(400 x 300 x 50 mm)
TG-HW-00095	TG-HW-000113	17" x 12" x 2"	(450 x 300 x 50 mm)
TG-HW-00096	TG-HW-000114	17" x 14" x 2"	(450 x 350 x 50 mm)
TG-HW-00097	TG-HW-000115	20 1/2" x 15 1/4" x 1 1/4"	(510 x 380 x 30 mm)
TG-HW-00098	TG-HW-000116	23 1/2" x 15 1/4" x 1 1/4"	(590 x 380 x 30 mm)
TG-HW-00099	TG-HW-000117	21" x 24" x 2"	(525 x 600 x 50 mm)

Mayo Tray

Length x Width x Height

TG-HW-000118	8" x 6" x 0.8"	(200 x 150 x 18 mm)
TG-HW-000119	8 1/2" x 6 1/4" x 0.8"	(215 x 160 x 18 mm)
TG-HW-000120	9 1/4" x 7 1/2" x 0.6"	(235 x 190 x 15 mm)
TG-HW-000121	11" x 7" x 0.4"	(280 x 180 x 10 mm)
TG-HW-000122	11 1/2" x 7 1/2" x 0.4"	(290 x 190 x 10 mm)
TG-HW-000123	12" x 8" x 0.9"	(300 x 200 x 20 mm)
TG-HW-000124	13 3/4" x 9 1/2" x 0.6"	(345 x 245 x 16 mm)
TG-HW-000125	14" x 9 1/2" x 0.4"	(350 x 240 x 10 mm)
TG-HW-000126	15 1/4" x 8 1/4" x 0.4"	(385 x 210 x 10 mm)
TG-HW-000127	15 1/4" x 10 1/2" x 0.6"	(385 x 265 x 16 mm)
TG-HW-000128	16" x 10 1/2" x 0.9"	(400 x 270 x 20 mm)
TG-HW-000130	16 3/4" x 11" x 0.4"	(420 x 280 x 10 mm)
TG-HW-000131	17" x 11 1/2" x 0.6"	(430 x 295 x 16 mm)
TG-HW-000132	19 1/4" x 13 3/4" x 0.6"	(485 x 320 x 16 mm)
TG-HW-000133	20 1/2" x 15 1/4" x 1 1/4"	(510 x 380 x 30 mm)
TG-HW-000134	23 1/2" x 15 1/4" x 1 1/4"	(590 x 380 x 30 mm)

Catheter Tray

Diameter x Height

TG-HW-000136	12" x 3" x 1 1/4"	(300 x 75 x 40 mm)
TG-HW-000137	14" x 3 1/2" x 1 1/2"	(350 x 90 x 50 mm)
TG-HW-000138	16" x 4" x 2"	(400 x 100 x 50 mm)
TG-HW-000139	18" x 4" x 2"	(450 x 100 x 50 mm)

Shallow Tray Grooved

Diameter x Height

TG-HW-000140	7 1/2" x 4 x 0.3"	(190 x 100 x 13 mm)
TG-HW-000141	8" x 4" x 0.3"	(200 x 100 x 13 mm)



Instruments Trays and Boxes

Trays Perforated

	Length x Diameter x Height	
TG-HW-000142	10" x 6" x 2"	(250 x 150 x 50 mm)
TG-HW-000143	12" x 6" x 2"	(300 x 150 x 50 mm)
TG-HW-000144	9" x 9" x 2"	(230 x 230 x 50 mm)
TG-HW-000145	19" x 10¼" x 2"	(480 x 260 x 50 mm)
TG-HW-000146	14" x 12" x 2"	(350 x 300 x 50 mm)



Rectangular Container Perforated

	Length x Diameter x Height	
TG-HW-000147	8½" x 4¾" x 2"	(220 x 120 x 50 mm)
TG-HW-000148	9" x 5" x 2"	(230 x 130 x 50 mm)
TG-HW-000149	9" x 6" x 2"	(230 x 150 x 50 mm)
TG-HW-000150	10" x 6" x 2"	(250 x 150 x 50 mm)
TG-HW-000151	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000152	12¾" x 3" x 2"	(325 x 75 x 50 mm)
TG-HW-000153	24" x 12" x 6"	(600 x 300 x 150 mm)



Instrument Box

	Diameter x Height	
TG-HW-000154	6¾" x 2¾" x 1¼"	(170 x 70 x 30 mm)
TG-HW-000155	7" x 3¼" x 1¾"	(180 x 80 x 40 mm)
TG-HW-000156	7" x 3½" x 1¾"	(180 x 90 x 40 mm)
TG-HW-000157	8" x 4" x 1¾"	(200 x 100 x 40 mm)
TG-HW-000158	8" x 4" x 2¾"	(200 x 100 x 60 mm)
TG-HW-000159	8½" x 4¾" x 2¾"	(220 x 120 x 60 mm)
TG-HW-000160	10" x 8" x 2¾"	(250 x 120 x 60 mm)
TG-HW-000161	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000162	12¾" x 6" x 4"	(320 x 150 x 60 mm)
TG-HW-000163	16" x 8" x 4"	(400 x 200 x 100 mm)
TG-HW-000164	20" x 8" x 4"	(500 x 200 x 100 mm)
TG-HW-000165	20" x 8¼" x 4¾"	(500 x 210 x 120 mm)



Instruments Box Perforated

	Diameter x Height	
TG-HW-000166	6¾" x 2¾" x 1¼"	(170 x 70 x 30 mm)
TG-HW-000167	7" x 3¼" x 1¾"	(180 x 80 x 40 mm)
TG-HW-000168	7" x 3½" x 1¾"	(180 x 90 x 40 mm)
TG-HW-000169	8" x 4" x 1¾"	(200 x 100 x 40 mm)
TG-HW-000170	8" x 4" x 2¾"	(200 x 100 x 60 mm)
TG-HW-000171	8½" x 4¾" x 2¾"	(220 x 120 x 60 mm)
TG-HW-000172	10" x 8" x 2¾"	(250 x 120 x 60 mm)
TG-HW-000173	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000174	12¾" x 6" x 4"	(320 x 150 x 60 mm)
TG-HW-000175	16" x 8" x 4"	(400 x 200 x 100 mm)
TG-HW-000176	20" x 8" x 4"	(500 x 200 x 100 mm)
TG-HW-000177	20" x 8¼" x 4¾"	(500 x 210 x 120 mm)



PRF Box, Endo Box Dental Cassette

Food Tray

Length x Diameter x Height

TG-HW-000178 16" x 12" x 1" (405 x 305 x 25 mm)



304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel

PRF System Box

Length x Diameter x Height

TG-HW-000179 10.47 x 4.09 x 2.24



304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel

Endodontic Box

Diameter x Height

TG-HW-000180 Blue 4 x 2.13 x 1.88 inches (10.2 x 5.4 x 4.8 cm)
 TG-HW-000181 Red 4 x 2.13 x 1.88 inches (10.2 x 5.4 x 4.8 cm)
 TG-HW-000182 Silver 4 x 2.13 x 1.88 inches (10.2 x 5.4 x 4.8 cm)

Dental bur block with 72 holes (36 holes for FG bur, point for High Speed, 36 holes for endodontic files, for K-Files or H-Files)
 Autoclavable under the high temperature of 135 degree and pressure of 0.22Mpa
 Material: Aluminum




Detachable Cassette

Lengthx width x Height

TG-HW-000183 Silver 14"x 9.5" x 1.4"
 TG-HW-000184 Black 14"x 9.5" x 1.4"

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel




Mesh Trays & Cassette

Mesh Cassette

Length x Diameter x Height

TG-HW-000185	(for 5 Instruments, Size 200 x 85 x 40 mm)
TG-HW-000186	(for 7 Instruments, Size Größe 200 x 105 x 40 mm)
TG-HW-000187	(for 10 Instruments, Size 200 x 145 x 35 mm)
TG-HW-000188	(for 15 Instruments, Size 200 x 200 x 35 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel




Instruments Cassette

Length x Diameter x Height

TG-HW-000189	(for 5 Instruments, Size 180 x 65 x 22 mm)
TG-HW-000190	(for 7 Instruments, Size Größe 180 x 90 x 22 mm)
TG-HW-000191	(for 10 Instruments, Size 185 x 135 x 22 mm)
TG-HW-000192	5 / 7 / 10 Instruments tray are same patren

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Mesh Trays W/Lid

Length x Width x Height

TG-HW-000193	40 x 40 X 20mm
TG-HW-000194	80 x 40 X 20mm
TG-HW-000195	160 x 40 X 20mm
TG-HW-000196	70 x 70 X 30mm
TG-HW-000197	80 x 80 X 40mm
TG-HW-000198	90 x 90 X 50mm
TG-HW-000199	180 x 130 x 20mm
TG-HW-000200	225 x 110 x 50mm
TG-HW-000201	270 x 135 x 40mm
TG-HW-000202	400 x 250 x 50mm

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel




Micro Mesh Tray

Length x Width x Height

TG-HW-000203	250 x 240 x 50mm
TG-HW-000204	250 x 250 x 50mm
TG-HW-000205	250 x 240 x 70mm
TG-HW-000206	250 x 250 x 50mm
TG-HW-000207	260 x 190 x 50mm
TG-HW-000208	380 x 250 x 50mm
TG-HW-000209	440 x 110 x 65mm
TG-HW-000210	480 x 250 x 50mm
TG-HW-000211	480 x 250 x 70mm
TG-HW-000212	500 x 150 x 100mm

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Surgical Container Boxes

Large/Small/Mini Micro Fragment Box

Length x Diameter x Height

TG-HW-000213 14.5"*8.5"*4.25"

Aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 2 sterilization trays . These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



DHS Box (Medium)

Length x Diameter x Height

TG-HW-000214 14.5"*8.5"*4.25"

Aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 3 sterilization trays . These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



DHS Box (Small)/Broken Set Box

Length x Width x Height

TG-HW-000215 14.5"*8.5"*3

Aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 2 sterilization trays . These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



ACL/PCL Box

Length x Width x Height

TG-HW-000216 19"*11"*4.25"

Aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 3 sterilization trays . These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



Surgical Container Boxes

Spine Fixation Box

Length x Diameter x Height

TG-HW-000217

21"*10"*5.2"

Anodized aluminum

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 3 sterilization trays . These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



Small Box With Tray

Length x Diameter x Height

TG-HW-000218

12"*8.5"*3"

Anodized aluminum

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 2 sterilization trays These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



Orthopedic Curette Box

Length x Width x Height

TG-HW-000219

11"*6"*2"

Anodized aluminum

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of gold box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 1 sterilization trays These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



90 Holes Pocket Size Golden Screw Box For 2mm

Length x Width x Height

TG-HW-000220

5.75"*3.75"*2.3"

Anodized aluminum

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of gold box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 2 sterilization trays These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



Surgical Container Boxes

Box For Holding Crocodile Forceps

Length x Diameter x Height

TG-HW-000221 12"*8"*4.5"

Aluminium

Sterilization screw box, aluminium box is used to carry crocodile forceps. Main body consist of box, silicone hoding & lid, trays and locks are produced from 1050 quality H24 aluminium box. These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



DHS DCS Box With 6.5,4.5mm Screw Rack

Length x Diameter x Height

TG-HW-000222 20.25"*9.5"*4.5

Anodized aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 3 sterilization trays
These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



Austin Moore Box

Length x Width x Height

TG-HW-000223 17"*8.5"*3"

Anodized aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, blue lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 2 sterilization trays .
These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.



Bookwalter Set Box

Length x Width x Height

TG-HW-000224 20*14*6

Anodized aluminium

Sterilization Container, aluminium box is used to carry and store the surgical tools and equipments. Main body consist of box, lid, trays and locks are produced from 1050 quality H24 aluminium box. It has 4 sterilization trays
These perforated orthopaedic surgical instruments sterilization boxes are pre-designed according user needs and before surgical operation it has to be sterilized covered with thick surgical textile and marked with indicators or labels.




Surgical Box, Microsurgery Box, Mats

Instrument Sterilization Box With Full Silicone Finger Mat

Width x Length x Height

TG-HW-000225	6 ³ / ₄ " x 2 ³ / ₄ " x 1 ¹ / ₄ "	(170 x 70 x 30 mm)
TG-HW-000226	7" x 3 ¹ / ₄ " x 1 ³ / ₄ "	(180 x 80 x 40 mm)
TG-HW-000227	7" x 3 ¹ / ₂ " x 1 ³ / ₄ "	(180 x 90 x 40 mm)
TG-HW-000228	8" x 4" x 1 ³ / ₄ "	(200 x 100 x 40 mm)
TG-HW-000229	8" x 4" x 2 ³ / ₄ "	(200 x 100 x 60 mm)
TG-HW-000230	8 ¹ / ₂ " x 4 ³ / ₄ " x 2 ³ / ₄ "	(220 x 120 x 60 mm)
TG-HW-000231	10" x 8" x 2 ³ / ₄ "	(250 x 120 x 60 mm)
TG-HW-000232	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000233	12 ³ / ₄ " x 6" x 4"	(320 x 150 x 60 mm)
TG-HW-000234	16" x 8" x 4"	(400 x 200 x 100 mm)
TG-HW-000235	20" x 8" x 4"	(500 x 200 x 100 mm)
TG-HW-000236	20" x 8 ¹ / ₄ " x 4 ³ / ₄ "	(500 x 210 x 120 mm)



304 Stainless steel high-quality steel
Autoclave/Ethylene Oxide/Dry heat 320°F (160°C)/Cold solution

Width x Length x Height

TG-HW-000237	10 1/2" x 6" x 1 1/8"
--------------	-----------------------



304 Stainless steel high-quality steel
Autoclave/Ethylene Oxide/Dry heat 320°F (160°C)/Cold solution

Microsurgery Instruments Box

Width x Length x Height

TG-HW-000238	10"x6"x1.5"
--------------	-------------



304 Stainless steel high-quality steel
Autoclave/Ethylene Oxide/Dry heat 320°F (160°C)/Cold solution

Finger Silicone Mat


Width x Length

TG-HW-000239	6 ³ / ₄ " x 2 ³ / ₄ "
TG-HW-000240	7" x 3 ¹ / ₄ "
TG-HW-000241	7" x 3 ¹ / ₂ "
TG-HW-000242	8" x 4" x
TG-HW-000243	8" x 4" x
TG-HW-000244	8 ¹ / ₂ " x 4 ³ / ₄ "
TG-HW-000245	10" x 8"
TG-HW-000246	12" x 8"
TG-HW-000247	12 ³ / ₄ " x 6"
TG-HW-000248	16" x 8"
TG-HW-000249	20" x 8"
TG-HW-000250	20" x 8 ¹ / ₄ "

Magnetic Silicone MAT

Width x Length

TG-HW-000251	12"x 15"
--------------	----------



Surgical Box, Microsurgery Box, Mats

Instruments Sterilization Tray Box

Width x Length x Height

TG-HW-000252	6 ³ / ₄ " x 2 ³ / ₄ " x 1 ¹ / ₄ "	(170 x 70 x 30 mm)
TG-HW-000253	7" x 3 ¹ / ₄ " x 1 ³ / ₄ "	(180 x 80 x 40 mm)
TG-HW-000254	7" x 3 ¹ / ₂ " x 1 ³ / ₄ "	(180 x 90 x 40 mm)
TG-HW-000255	8" x 4" x 1 ³ / ₄ "	(200 x 100 x 40 mm)
TG-HW-000256	8" x 4" x 2 ³ / ₄ "	(200 x 100 x 60 mm)
TG-HW-000257	8 ¹ / ₂ " x 4 ³ / ₄ " x 2 ³ / ₄ "	(220 x 120 x 60 mm)
TG-HW-000258	10" x 8" x 2 ³ / ₄ "	(250 x 120 x 60 mm)
TG-HW-000259	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000260	12 ³ / ₄ " x 6" x 4"	(320 x 150 x 60 mm)
TG-HW-000261	16" x 8" x 4"	(400 x 200 x 100 mm)
TG-HW-000262	20" x 8" x 4"	(500 x 200 x 100 mm)
TG-HW-000263	20" x 8 ¹ / ₄ " x 4 ³ / ₄ "	(500 x 210 x 120 mm)

Autoclavable Sterile Tray Box used to hold Orthopedic, Microsurgery, Cardiac surgery Ophthalmic Surgery Instruments Available in customize sizes and colors .

Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize

Instruments Sterilization Container Box with Silicone Mat

Width x Length x Height

TG-HW-000264	6 ³ / ₄ " x 2 ³ / ₄ " x 1 ¹ / ₄ "	(170 x 70 x 30 mm)
TG-HW-000265	7" x 3 ¹ / ₄ " x 1 ³ / ₄ "	(180 x 80 x 40 mm)
TG-HW-000266	7" x 3 ¹ / ₂ " x 1 ³ / ₄ "	(180 x 90 x 40 mm)
TG-HW-000267	8" x 4" x 1 ³ / ₄ "	(200 x 100 x 40 mm)
TG-HW-000268	8" x 4" x 2 ³ / ₄ "	(200 x 100 x 60 mm)
TG-HW-000269	8 ¹ / ₂ " x 4 ³ / ₄ " x 2 ³ / ₄ "	(220 x 120 x 60 mm)
TG-HW-000270	10" x 8" x 2 ³ / ₄ "	(250 x 120 x 60 mm)
TG-HW-000271	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000272	12 ³ / ₄ " x 6" x 4"	(320 x 150 x 60 mm)
TG-HW-000273	16" x 8" x 4"	(400 x 200 x 100 mm)
TG-HW-000274	20" x 8" x 4"	(500 x 200 x 100 mm)
TG-HW-000275	20" x 8 ¹ / ₄ " x 4 ³ / ₄ "	(500 x 210 x 120 mm)

Autoclavable Sterile Tray Box used to hold Orthopedic, Microsurgery, Cardiac surgery Ophthalmic Surgery Instruments Available in customize sizes and colors .

Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize

Instruments Sterilization Container Box

Width x Length x Height

TG-HW-000276	6 ³ / ₄ " x 2 ³ / ₄ " x 1 ¹ / ₄ "	(170 x 70 x 30 mm)
TG-HW-000277	7" x 3 ¹ / ₄ " x 1 ³ / ₄ "	(180 x 80 x 40 mm)
TG-HW-000278	7" x 3 ¹ / ₂ " x 1 ³ / ₄ "	(180 x 90 x 40 mm)
TG-HW-000279	8" x 4" x 1 ³ / ₄ "	(200 x 100 x 40 mm)
TG-HW-000280	8" x 4" x 2 ³ / ₄ "	(200 x 100 x 60 mm)
TG-HW-000281	8 ¹ / ₂ " x 4 ³ / ₄ " x 2 ³ / ₄ "	(220 x 120 x 60 mm)
TG-HW-000282	10" x 8" x 2 ³ / ₄ "	(250 x 120 x 60 mm)
TG-HW-000283	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000284	12 ³ / ₄ " x 6" x 4"	(320 x 150 x 60 mm)
TG-HW-000285	16" x 8" x 4"	(400 x 200 x 100 mm)
TG-HW-000286	20" x 8" x 4"	(500 x 200 x 100 mm)
TG-HW-000287	20" x 8 ¹ / ₄ " x 4 ³ / ₄ "	(500 x 210 x 120 mm)

Autoclavable Sterile Tray Box used to hold Orthopedic, Microsurgery, Cardiac surgery Ophthalmic Surgery Instruments Available in customize sizes and colors .

Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize


Instruments Sterilization Container Box

Width x Length x Height

TG-HW-000288	6 ³ / ₄ " x 2 ³ / ₄ " x 1 ¹ / ₄ "	(170 x 70 x 30 mm)
TG-HW-000289	7" x 3 ¹ / ₄ " x 1 ³ / ₄ "	(180 x 80 x 40 mm)
TG-HW-000290	7" x 3 ¹ / ₂ " x 1 ³ / ₄ "	(180 x 90 x 40 mm)
TG-HW-000291	8" x 4" x 1 ³ / ₄ "	(200 x 100 x 40 mm)
TG-HW-000292	8" x 4" x 2 ³ / ₄ "	(200 x 100 x 60 mm)
TG-HW-000293	8 ¹ / ₂ " x 4 ³ / ₄ " x 2 ³ / ₄ "	(220 x 120 x 60 mm)
TG-HW-000294	10" x 8" x 2 ³ / ₄ "	(250 x 120 x 60 mm)
TG-HW-000295	12" x 8" x 2"	(300 x 200 x 50 mm)
TG-HW-000296	12 ³ / ₄ " x 6" x 4"	(320 x 150 x 60 mm)

Autoclavable Sterile Tray Box used to hold Orthopedic, Microsurgery, Cardiac surgery Ophthalmic Surgery Instruments Available in customize sizes and colors .

Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize



Surgical Box, Microsurgery Box, Mats

Box For Holding Crocodile Forceps

Width x Length x Height

TG-HW-000297 12"X8"X4.5" (300mmX200mmX109mm)

Basket case ultrasonic Cleaning Washing Trays Sterilization by all standard methods Autoclave
Available in customize sizes and colors .
Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize




Ophthalmic Sterilization Tray Box with Silicone Mat

Width x Length x Height

TG-HW-000298 (150mmX120mmX28mm)
TG-HW-000299 (195mmX104mmX27mm)
TG-HW-000300 (260mmX152mmX28mm)

Basket case ultrasonic Cleaning Washing Sterilization Trays use for Microsurgery,
Ophthalmic Instruments Sterilization by all standard methods Autoclave
Available in customize sizes and colors .
Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize




Aluminum Sterilization Box Double Tray With Mat

Width x Length x Height

TG-HW-000301 12³/₄" x 6" x 4" (320 x 150 x 60 mm)
TG-HW-000302 16" x 8" x 4" (400 x 200 x 100 mm)
TG-HW-000303 20" x 8" x 4" (500 x 200 x 100 mm)
TG-HW-000304 20" x 8¹/₄" x 4³/₄" (500 x 210 x 120 mm)

Autoclavable Sterile Tray Box used to hold Orthopedic, Microsurgery, Cardiac surgery
Ophthalmic Surgery Instruments Available in customize sizes and colors .
Produced from 1050 Quality H24 Aluminium . 130° to 135° Temperature Sterilize




Standard Tray, Deep base, Insert Tray, Lid & 2 Silicone Mats

Width x Length x Height

TG-HW-000305 6.0 X 10.0 X 1.5"

These features make PST Plastic Sterilization Trays better for sterilizing your surgical instruments.
Flexibility of layout. Either finger-tip mats or special, "plug-in" inserts for individual instruments.
Unique, "grid" system in all bases, trays and lids, making it easy to install the mat or insert required to protect your instruments properly.
Stackable trays. Integral feet in bases and corner ridges on lids allow for easier handling in autoclaves and storage.
All bases have pebble-type surface to keep mats from sticking.
Dimensions of each tray are designed to accommodate a variety of instrument sizes and shapes.
Sterilizable by all standard methods of sterilization.



Needle Cases, Petri Dishes, Cotton Jars, Spirit Lamp


Needle Cases

Complete with perforated interior box

	Diameter x Height
TG-HW-000306	Ø 65 mm x 15 mm

Interior box perforated lid and bottom

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Petri Dishes

	Diameter x Height
TG-HW-000307	2¾" x 0.9" (60 x 20 mm)
TG-HW-000308	3½" x 0.6" (90 x 15 mm)
TG-HW-000309	4" x 0.9" (100 x 20 mm)
TG-HW-000310	4¾" x 0.9" (120 x 20 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Spirit Lamp

	Size
TG-HW-000311	8cm

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Cotton and Catgut Jars

	Diameter x Height
TG-HW-000312	3" x 3" (75 x 75 mm)
TG-HW-000313	4" x 4" (100 x 100 mm)
TG-HW-000314	5" x 5" (125 x 125 mm)
TG-HW-000315	6" x 6" (150 x 150 mm)
TG-HW-000316	7" x 7" (180 x 180 mm)
TG-HW-000317	8" x 8" (200 x 200 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Beakers

	Diameter x Height	Liter
TG-HW-000318	3" x 3" (75 x 75 mm)	0.25
TG-HW-000319	4" x 3½" (100 x 85 mm)	0.50
TG-HW-000320	5" x 4½" (115 x 100 mm)	0.75
TG-HW-000321	5" x 5" (125 x 125 mm)	1.00

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Dressing & Cotton Drums

Dressing Drums Round, spot welded

	Diameter x Height	
TG-HW-000321	3¼" x 3¼"	(80 x 80 mm)
TG-HW-000322	4¾" x 4¾"	(120 x 120 mm)
TG-HW-000323	5" x 5"	(125 x 125 mm)
TG-HW-000324	5½" x 5½"	(140 x 140 mm)
TG-HW-000325	6" x 6"	(150 x 150 mm)
TG-HW-000326	6¼" x 4¾"	(160 x 120 mm)
TG-HW-000327	6¼" x 6"	(160 x 150 mm)
TG-HW-000328	6¼" x 6¼"	(160 x 160 mm)
TG-HW-000329	7½" x 6¾"	(190 x 170 mm)
TG-HW-000330	7½" x 7½"	(190 x 190 mm)
TG-HW-000331	8" x 8"	(200 x 200 mm)
TG-HW-000332	9" x 9"	(230 x 230 mm)
TG-HW-000333	9" x 13¼"	(230 x 330 mm)
TG-HW-000334	9½" x 6¼"	(240 x 160 mm)
TG-HW-000335	9½" x 9½"	(240 x 240 mm)
TG-HW-000336	9½" x 11½"	(240 x 290 mm)
TG-HW-000337	9½" x 13½"	(240 x 340 mm)
TG-HW-000338	10" x 10"	(250 x 250 mm)
TG-HW-000339	11" x 11"	(280 x 280 mm)
TG-HW-000340	11½" x 6¼"	(290 x 160 mm)
TG-HW-000341	11¼" x 9½"	(290 x 240 mm)
TG-HW-000342	12" x 12"	(300 x 300 mm)
TG-HW-000343	13½" x 6¼"	(340 x 160 mm)
TG-HW-000344	13½" x 9½"	(340 x 240 mm)
TG-HW-000345	15¼" x 12"	(380 x 300 mm)
TG-HW-000346	15¼" x 15¼"	(380 x 380 mm)
TG-HW-000347	15½" x 7½"	(390 x 190 mm)
TG-HW-000348	15½" x 11½"	(390 x 290 mm)



304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel

Dressing Drum Jointless

	Diameter x Height	
TG-HW-000349	7½" x 5½"	(190 x 140 mm)
TG-HW-000350	9½" x 4¾"	(240 x 120 mm)
TG-HW-000351	9½" x 6¼"	(240 x 160 mm)
TG-HW-000352	13½" x 6¼"	(340 x 160 mm)
TG-HW-000353	15½" x 5½"	(390 x 140 mm)
TG-HW-000354	15½" x 6¼"	(390 x 160 mm)
TG-HW-000355	13½" x 9"	(340 x 230 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Jugs, Feeding & Septum Cup, Strainer

Graduated Jugs

	Liter
TG-HW-000356	0.25
TG-HW-000357	0.50
TG-HW-000358	1.00
TG-HW-000359	2.00

High-quality stainless steel measuring cups with practical measuring scale in timeless design.

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Feeding Cup

	Diameter x Height
TG-HW-000360	4" x 2" (100 x 50 mm)

Feeding Cup, diameter 4" (100 mm), height 2" (50 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Funnels

	Width x Length x Height
TG-HW-000361	4" x 4" (100 x 100 mm)
TG-HW-000362	5" x 5" (150 x 150 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Ball Tea Strainer Fine Mesh Stainless Steel, With Spring Closure

	Ring Sizes	Handle
TG-HW-000363	5 CM INNER DIAMETER	10cm
TG-HW-000364	7.5 CM INNER DIAMETER	10cm
TG-HW-000365	9.5 CM INNER DIAMETER	10cm

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Bed Pans and Urinal Pot

Bed Pan

	Length x Width x Height	for
TG-HW-000366	15¼" x 9" x 3¾" (380 x 230 x 95 mm)	Child
TG-HW-000367	14" x 10½" x 4" (350 x 270 x 100 mm)	Adult

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Female Urinal Jar

	Length x Width x Height
TG-HW-000368	9" x 6" x 3½" (230 x 150 x 85 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Male Urinal Jar

	Length x Width x Height
TG-HW-000369	9" x 6" x 3½" (230 x 150 x 85 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Bed Pans

For Fracture Patients

With cover	Length x Width x Height
TG-HW-000370	9" x 6" x 3½" (230 x 150 x 85 mm)

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel



Hospital Bucket & Urinal, Can

Bed Pan Urinal Jug type

Diameter X Height

TG-HW-000371 4½" x 9" (85 x 230 mm)



304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel

Female Urinal Jar

Diameter x Height

Liter

TG-HW-000372	4" x 5"	(100 x 125 mm)	1.0
TG-HW-000373	5" x 6"	(125 x 150 mm)	1.5
TG-HW-000374	5" x 7"	(125 x 175 mm)	2.0



304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel

Buckets

Without cover

With cover

Liter

TG-HW-000375	TG-HW-00082	10.0
TG-HW-000376	TG-HW-00082	13.0
TG-HW-000377	TG-HW-00082	18.0
TG-HW-000378	TG-HW-00082	24.0



304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel

Kick bucket and stand

With cover

With cover

Liter

TG-HW-000379	TG-HW-00082	10.0
TG-HW-000380	TG-HW-00082	13.0
TG-HW-000381	TG-HW-00082	18.0
TG-HW-000382	TG-HW-00082	24.0

304 Stainless steel is a strong material and stronger than others. It is easy to assemble and high-quality steel





Tahir Group Of Industries

HOSPITAL & MEDICAL HOLLOWWARE MANUFACTURER



ISO



Medicines & Healthcare products Regulatory Agency

