



Sections

- 01 Single Use Laryngoscope
- 02 Flexible Laryngoscope
- 03 Reusable Laryngoscope
- 04 Rechargeable ENT Diagnostic Set
- 05 Mini Fiber Optic Oscope
- 06 Metal Oscope Ophthalmoscope

Meeting The New Era Challenges
By Enabling Customer Satisfaction





01

Single Use
Laryngoscope

SINGLE USE LARYNGOSCOPE



Flexible
Blade Fiber Optic

LED BLUB

TGIND:5142

High Grade Surgical Steel
Spatula Provides With Unbreakable Plastic Base
Unique Heel Design Prevents Blade From Touching The Handle
Clearly Marked & Packed Single Patient Use
Fiber Optic Tube Black Jacket End Lit Fiber
Smooth Edges, Rounded Tip & High Quality Surface Finsih
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5142	1	76
TGIND:5142	2	88
TGIND:5142	3	112
TGIND:5142	4	130



Flexible
Blade Fiber Optic

LED BLUB

TGIND:5143

High Grade Surgical Steel
Spatula Provides With Unbreakable Plastic Base
Unique Heeal Design Prevents Blade From Touching The Handle
Clearly Marked & Packed Single Patient Use
Smooth Edges, Rounded Tip & High Quality Surface Finsih
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5143	1	76
TGIND:5143	2	88
TGIND:5143	3	112
TGIND:5143	4	130

SINGLE USE LARYNGOSCOPE



TGIND:5144

High Grade Surgical Steel
Spatula Provides With Unbreakable Plastic Base
Unique Heel Design Prevents Blade From Touching
The Handle Clearly Marked & Packed Single Patient Use
Smooth Edges, Rounded Tip & High Quality Surface Finsih
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5144	1	76
TGIND:5144	2	88
TGIND:5144	3	112
TGIND:5144	4	130



TGIND:5145

High Grade Surgical Steel
Spatula Provides With Unbreakable Plastic Base
Unique Heel Design Prevents Blade From Touching
The Handle Clearly Marked & Packed Single Patient Use
Smooth Edges, Rounded Tip & High Quality Surface Finsih
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5145	1	76
TGIND:5145	2	88
TGIND:5145	3	112
TGIND:5145	4	130

SINGLE USE LARYNGOSCOPE



TGIND:5146

High Grade Surgical Steel
Spatula Provides With Unbreakable Plastic Base
Unique Heel Design Prevents Blade From Touching
The Handle Clearly Marked & Packed Single Patient Use
Smooth Edges, Rounded Tip & High Quality Surface Finish
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5146	1	76
TGIND:5146	2	88
TGIND:5146	3	112
TGIND:5146	4	130



TGIND:5147

High Grade Surgical Steel
Spatula Provides With Unbreakable Plastic Base
Unique Heel Design Prevents Blade From Touching
The Handle Clearly Marked & Packed Single Patient Use
Smooth Edges, Rounded Tip & High Quality Surface Finish
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5147	1	76
TGIND:5147	2	88
TGIND:5147	3	112
TGIND:5147	4	130

DISPOSABLE LARYNGOSCOPE



TGIND:5148

Fiber Optic Disposable Laryngoscope Plastic Handle with Jointed Aluminium Head
Blades Miller 00, 0, 1, 2, 3,
Packing PU Leather Pouch



TGIND:5149

Fiber Optic Disposable Laryngoscope Plastic Handle with Jointed Aluminium Head
Plastic Blades Macintosh 0, 1, 2, 3
Packing PU Leather Pouch

DISPOSABLE LARYNGOSCOPE



TGIND:5150

TGIND:5151

TGIND:5152

TGIND:5153

TGIND:5154


Fiber Optic Disposable Laryngoscope Aluminium Handle with
Jointed Plastic Head

Macintosh **Blade Size: 0, 1, 2, 3, 4**

Packing PU Leather Pouch



SINGLE USE LARYNGOSCOPE



TGIND:5155

Tip to Spatula 303/304 Surgical Stainless Steel Heel Silver
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearing Gives Secure Fit and Firm Support
Fiber Optic Tube Black Jacket End Lit Fiber
Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5155	1	76
TGIND:5155	2	88
TGIND:5155	3	112
TGIND:5155	4	130




TGIND:5156

Tip to Spatula 303/304 Surgical Stainless Steel Heel Silver
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearing Gives Secure Fit and Firm Support
Fiber Optic Acrylic Rod
Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5156	00	45
TGIND:5156	0	55
TGIND:5156	1	80
TGIND:5156	2	130
TGIND:5156	3	170
TGIND:5156	4	180


SINGLE USE LARYNGOSCOPE



TGIND:5157

Tip to Spatula 303/304 Surgical Stainless Steel Heel Plastic (ABS) Smooth Edges, Rounded Tip & High Quality Surface Finish Two Spring Ball Bearing Gives Secure Fit and Firm Support Fiber Optic Acrylic Rod Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5157	1	76
TGIND:5157	2	88
TGIND:5157	3	112
TGIND:5157	4	130



TGIND:5158

Spatula 303/304 Surgical Stainless Steel Heel Plastic (ABS) Smooth Edges, High Quality Surface Finish Two Spring Ball Bearing Gives Secure Fit and Firm Support Fiber Optic Acrylic Rod Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5158	00	45
TGIND:5158	0	55
TGIND:5158	1	80
TGIND:5158	2	130
TGIND:5158	3	170
TGIND:5158	4	180

SINGLE USE LARYNGOSCOPE



TGIND:5159

Tip to Spatula 303/304 Surgical Stainless Steel Heel Silver
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearing Gives Secure Fit and Firm Support
Fiber Optic Tube Black Jacket End Lit Fiber
Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5159	1	76
TGIND:5159	2	88
TGIND:5159	3	112
TGIND:5159	4	130



TGIND:5160

Tip to Spatula 303/304 Surgical Stainless Steel Heel Silver
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearing Gives Secure Fit and Firm Support
Fiber Optic Acrylic Rod
Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5160	1	76
TGIND:5160	2	88
TGIND:5160	3	112
TGIND:5160	4	130

SINGLE USE LARYNGOSCOPE



TGIND:5161

Tip to Heel 303/304 Surgical Stainless Steel
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearing Gives Secure Fit and Firm Support
Fiber Optic Tube Black Jacket End Lit Fiber
Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5161	1	82
TGIND:5161	2	92
TGIND:5161	3	116
TGIND:5161	4	135




TGIND:5162

Tip to Heel 303/304 Surgical Stainless Steel
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearing Gives Secure Fit and Firm Support
With Plastic Sleeve 2.5v LED 5000 Lux Brightness
Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5162	1	82
TGIND:5162	2	92
TGIND:5162	3	116
TGIND:5162	4	135

SINGLE USE LARYNGOSCOPE



TGIND:5163

Tip to Spatula 303/304 Surgical Stainless Steel Unique Heel Design Smooth Edges, Rounded Tip & High Quality Surface Finish Two Spring Ball Bearing Gives Secure Fit and Fim Support 2.5 Single & Double LED Design Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5163	1	83
TGIND:5163	2	95
TGIND:5163	3	120
TGIND:5163	4	138



TGIND:5164

Tip to Spatula 303/304 Surgical Stainless Steel Unique Heel Design Smooth Edges, High Quality Surface Finish Two Spring Ball Bearing Gives Secure Fit and Fim Support 2.5 Single & Double LED Design Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5164	00	45
TGIND:5164	0	55
TGIND:5164	1	82
TGIND:5164	2	130
TGIND:5164	3	175
TGIND:5164	4	185

REUSABLE LARYNGOSCOPE

Macintosh F/O Fix Tube (GS 100)



TGIND:5165

303/304 Surgical Stainless Steel
 Smooth Edges, Rounded Tip & High Quality Surface Finish
 Provides with 3mm Fiber Optic Offers No Black Marks
 Fully Autoclavable & Easier to Clean
 Satin Finish

Macintosh	Fig	Working Length(mm)
TGIND:5165	1	80
TGIND:5165	2	92
TGIND:5165	3	116
TGIND:5165	4	135

SINGLE USE LARYNGOSCOPE

Macintosh	Fig	Working Length(mm)
TGIND.5166	1	80
TGIND.5166	2	92
TGIND.5166	3	116
TGIND.5166	4	133

TGIND:5166

Tip to Spatula 303/304 Surgical Stainless Steel
 Unique Heel Design Smooth Edges, Rounded Tip
 & High Quality Surface Finish Two Spring Ball
 Bearing Gives Secure Fit and Firm Support
 2.5 Single LED Design 5000 Lux Brightness
 Satin Finish



Conventional Standard



Medium

Small

Stubby

Fiber Standard



Medium

Small

Stubby

Single Use Handle




Conventional

Conventional


Fiber Optic

SINGLE USE LARYNGOSCOPE WTH BOX


Single Use Laryngoscope Set
With 4 Macintosh Blades 1, 2, 3, 4




Single Use Laryngoscope Set
With 4 Macintosh Blades 1, 2, 3, 4



Single Use Laryngoscope Set
With 3 Flexiable Macintosh Blades 1, 2, 3, 4



Single Use Laryngoscope Set
With 4 Miller Blades 00, 0, 1, 2, 3, 4, 6, 7





Flexible Fiber Optic Laryngoscope

02



REUSABLE LARYNGOSCOPE



Macintosh Flexible
Fix Tube Fiber Optic
Cold Light Guide

LED/HALOGEN BLUB

TGIND:5171

High Grade Surgical Steel
Smooth Edges, Rounde Tip & High Quality Surface Finish
Provides With 3mm Fiber Opitic Offer More No Black Marks
Fully Autoclavable & Easier to Clean
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5171	1	70
TGIND:5171	2	88
TGIND:5171	3	112
TGIND:5171	4	134



Macintosh Flexible
Fix Tube Fiber Optic
Cold Light Guide

LED/HALOGEN BLUB

TGIND:5172

High Grade Surgical Steel
Smooth Edges, Rounde Tip & High Quality Surface Finish
Three Spring Ball Bearings Gives Secure Fit & Support
Supplied With LED
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5172	1	92
TGIND:5172	2	113
TGIND:5172	3	135
TGIND:5172	4	15

REUSABLE LARYNGOSCOPE



Macintosh Flexible
Fix Tube Fiber Optic

LED/HALOGEN BLUB

TGIND:5173

High Grade Surgical Steel
Smooth Edges, Rounde Tip & High Quality Surface Finish
Large Fiber Bundles Offer More Light No Black Marks
Fully Autoclavable & Easier to Clean
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5173	1	72
TGIND:5173	2	88
TGIND:5173	3	112
TGIND:5173	4	134



Macintosh Flexible
Fix Tube Fiber Optic

LED/HALOGEN BLUB

TGIND:5174

High Grade Surgical Steel
Smooth Edges, Rounde Tip & High Quality Surface Finish
Removeable Light Guides Design Allow Easier Cleaning &
Extends Blade Life Provide Wih 4mm Fiber Optic Bundles
for Light Pass Gives True Tissue Color
Fully Autoclavable & Easier to Clean
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5174	1	80
TGIND:5174	2	92
TGIND:5174	3	116
TGIND:5174	4	135



03

Reusable
Laryngoscope

REUSABLE LARYNGOSCOPE

Macintosh Fix Tube (Fiber Optic)
Cold Light Guide



TGIND:5175

303/304 Surgical Stainless Steel
Smooth Edges, Rounded Tip & High Quality Surface Finish
Provide With 3mm Fiber Optic Bundles For Light Pass Gives
Truse Tissue Color
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5175	1	70
TGIND:5175	2	88
TGIND:5175	3	112
TGIND:5175	4	134

Miller Fix Tube (Fiber Optic)
Cold Light Guide



TGIND:5176

303/304 Surgical Stainless Steel
Smooth Edges, High Quality Surface Finish
Provide with 3mm Fiber Optic Bundles for Light Pass Gives
True Tissue Color
Fully Autoclavable & Easier to Clean
Satin Finish

Miller	Fig	Working Length(mm)
TGIND:5176	00	45
TGIND:5176	0	57
TGIND:5176	1	80
TGIND:5176	2	130
TGIND:5176	3	170
TGIND:5176	4	180

REUSABLE LARYNGOSCOPE

Macintosh Interchangable Blade
(Fiber Optic) Cold Light Guide



TGIND:5177

303/304 Surgical Stainless Steel
Smooth Edges, Rounde Tip & High Quality Surface Finish
Removable Light Guides Design Alllow Easier Cleaning &
Extends Blade Life Provide With 4mm Fiber Optic Bundles
For Light Pass Gives Truse Tissue Color
Stin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5177	1	80
TGIND:5177	2	92
TGIND:5177	3	116
TGIND:5177	4	135

Miller Interchangable Blade
(Fiber Optic) Cold Light Guide



TGIND:5178

303/304 Surgical Stainless Steel
Smooth Edges, High Quality Surface Finish
Provide with 4mm Fiber Optic Bundles for Light Pass Gives
True Tissue Color
Fully Autoclavable & Easier to Clean
Satin Finish

Miller	Fig	Working Length(mm)
TGIND:5178	00	45
TGIND:5178	0	55
TGIND:5178	1	80
TGIND:5178	2	130
TGIND:5178	3	170
TGIND:5178	4	180

REUSABLE LARYNGOSCOPE

Macintosh Conventional



TGIND:5179

303/304 Surgical Stainless Steel
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearings Gives Secure
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5179	0	60
TGIND:5179	1	68
TGIND:5179	2	92
TGIND:5179	3	113
TGIND:5179	4	135
TGIND:5179	5	155

Miller Conventional



TGIND:5180

303/304 Surgical Stainless Steel
Smooth Edges, High Quality Surface Finish
Three Spring Ball Bearings Gives Secue Fit & Support
Supplied With LED
Satin Finish

Miller	Fig	Working Length(mm)
TGIND:5180	00	60
TGIND:5180	0	70
TGIND:5180	1	92
TGIND:5180	2	113
TGIND:5180	3	135
TGIND:5180	4	155

REUSABLE LARYNGOSCOPE

Integrated Laryngoscope (Fiber Optic)
Cold Light Guide



TGIND:5181

303/304 Surgical Stainless Steel
Smooth Edges, Round Tip & High Quality Surface Finish
Large Fiber Bundle Offers More Light no Black Marks
Fully Autoclavable & Easier to Clean
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5181	1	72
TGIND:5181	2	90
TGIND:5181	3	112
TGIND:5181	4	135

Miller Fix Tube HBT (Fiber Optic)
Cold Light Guide



TGIND:5182

303/304 Surgical Stainless Steel
Smooth Edges, High Quality Surface Finish
Provide with 3mm Fiber Optic Offers No Black Marks
Fully Autoclavable & Easier to Clean
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5182	00	45
TGIND:5182	0	55
TGIND:5182	1	80
TGIND:5182	2	130
TGIND:5182	3	170
TGIND:5182	4	130

REUSABLE LARYNGOSCOPE

Macintosh Conventional Eco



TGIND:5183

303/304 Surgical Stainless Steel
Smooth Edges, Rounded Tip & High Quality Surface Finish
Three Spring Ball Bearings Gives Secure Fit & Firm Support Supplied
With LED 2.5V 5000 Lux
Satin Finish

MACINTOSH	Fig	Working Length(mm)
TGIND:5183	1	68
TGIND:5183	2	92
TGIND:5183	3	113
TGIND:5183	4	135

Miller Conventional Eco



TGIND:5184

303/304 Surgical Stainless Steel
Smooth Edges, High Quality Surface Finish
Three Spring Ball Bearings Gives Secure Fit & Firm Support Supplied
with LED 2.5V 5000 Lux
Satin Finish

Miller	Fig	Working Length(mm)
TGIND:5184	00	45
TGIND:5184	0	57
TGIND:5184	1	80
TGIND:5184	2	130
TGIND:5184	3	170
TGIND:5184	4	180

REUSABLE LARYNGOSCOPE



TGIND:5185

Integrated Fiber Optic Laryngoscope Set LED Light
With battery handle
4 Blades 1, 2, 3, 4 with Transparent Box



TGIND:5186

Macintosh Fiber Optic Laryngoscope Set LED Light
With battery handle
4 Blades 1, 2, 3, 4 with Transparent Box

REUSABLE LARYNGOSCOPE



TGIND:5187

Conventional
Miller & Macintosh Pediatric Laryngoscope Set
6 Blades With Transparent Box Red

Macintosh Blade 0,
Macintosh Blade 1,
Macintosh Blade 2
Miller Blade 0,
Miller Blade 1,
Miller Blade 2



TGIND:5188

Conventional
Miller Laryngoscope Set 6 straight Blades
With Transparent Box Red

Size 0, Premie; Overall Length: 80 mm.
Size 1, Infant; Overall Length: 105 mm.
Size 2, Child; Overall Length: 155 mm.
Size 3, Adult; Overall Length: 195 mm.
Size 4, Large Adult; Overall Length: 205 mm.
Size 5, Extra Large Adult; Overall Length: 220 mm.



04



Rechargeable
Diagnostic ENT Set

OPHTHALMOSCOPE OTOSCOPE RECHARGEABLE

TGIND:5189

Fiber Optic Otoscope Ophthalmoscope
Adult & Pediatric Rechargeable Dekstop



TGIND:5190


Fiber Optic Otoscope + Ophthalmoscope
Pediatric With Wall-Mounted Charger



HANDLES CHARGER

TGIND:5191

Fiber Optic Laryngoscope
Adult & Pediatric Handles Charger Dekstop



TGIND:5192

Rechargeable Laryngoscope Handles are designed to allow even simpler, faster and safer hygienic reprocessing of the handle has no complex disassembling and assembling of batteries and light unit. This saves the technical service personnel with a great deal of time, and reduces the training and maintenance effort. Handle comes with 3.5V Ni-MH 1000mAh rechargeable battery, can be recharge with Pro Charger.



OTOSCOPE RECHARGEABLE

TGIND:5193

Conventional Otoscope
Rechargeable LED light
3 Cannula With UBS Cable
& Pouch



TGIND:5194

Fiber Optic Mini Otoscope Rechargeable

5 Pcs of adult spatulas.
5 Pcs of Infant/Child spatulas.
1 Extra LED Bulb.
1 USB Cable
1 Pouch



OPHTHALMOSCOPE OTOSCOPE



TGIND:5210

ENT Ophthalmoscope Mini Fiber Optic Otoscope Set
for adult & pediatric

Otoscope Scope with 2 Light,
Ophthalmoscope
Pneumatic Air Pump for Insufflation Ear Infection Detector
8 Disposable Specula
packing in transparent box



TGIND:5211

Diagnostic ENT Ophthalmoscope Otoscope Set

Ophthalmoscope
Otoscope
Plastic Reusable Specula
2 Plastic Reusable Specula
Chrome Plated And Knurled Handle
2 Laryngeal mirrors
Charger
Rechargeable Battery

DIAGNOSTIC ENT SET

Diagnostic ENT Examination Set

TGIND:5208

- 1 Handle (takes 2 "C" batteries not included)
- 3 Ear Specula (DIFFERENT SIZES)
- 1 Bent Arm Illuminator
- 1 Nasal Adapter
- 1 Head Lens Wheel
- 2 Mirrors #4 #5
- 4 Free Extra Bulbs
- 1 Tongue Depressor
- 20 spear mini specula size 2.5 mm to 4.0 mm
- All Packed in a Zipper Pouch as Shown In the Picture is of the Actual Product



cc

Diagnostic ENT Set

TGIND:5209

ENT Diagnostic Set is illuminated by a 2.7 V vacuum lamp. The fluted handle surface is hygienic and allows for a secure grip. The ENT Set plug into wall plugs for quick and easy charging of the batteries. The Otoscope and Ophthalmoscope set includes: Otoscope head, ophthalmoscope head, expandable nasal speculum, bent arm illuminator, 2 laryngeal mirrors, tongue blade holder, replacement otoscope bulb, and hard protective case.



OPHTHALMOSCOPE OTOSCOPE

TGIND:5212

ENT Ophthalmoscope Otoscope Nasal Larynx Diagnostic Set

Otoscope:

Chrome plated Otoscope head with removable 4x magnifying lens and inufflator fitting

Autoclavable polypropylene ear specula in 2.5mm, 3.5 mm and 4.5 mm

14 speculas for children

Nasal speculum adaptor with manual screw

opening for instrument connected to the Otoscope head

Ophthalmoscope:

Ophthalmoscope head with lens wheel featuring 25 corrective lens from - 25 to + 40 diopters

Bent Arm Illuminator:

Chrome plated bent arm larynx with a bulb One tongue depressor adaptor Stainless Steel
2x 4.0 mm mirrors





05

Fiber Optic Mini Otoscope

FIBER OPTIC OTOSCOPE

TGIND:5196

Fiber Mini Otoscope Kit

- 1 Otoscope
- 1 Ophthalmoscope
- 5 Pcs of adult spatulas.
- 5 Pcs of Infant/Child spatulas.
- 2 Extra LED Bulb.
- 1 Pouch



TGIND:5197

Fiber Mini Otoscope Kit

- 1 Otoscope
- 1 Ophthalmoscope
- 5 Pcs of adult spatulas.
- 5 Pcs of Infant/Child spatulas.
- 2 Extra LED Bulb.
- 1 Pouch



MINI FIBER OPTIC OTOSCOPE

TGIND:5195

Fiber Optical
(Pocket Style)
2.5v LED

Modern pocket otoscope compact design. Handy and powerful. Distal Fiber-Optic (F.O.) illumination 4mm bundle with Enhanced LED Technology. Attachment clip with integrated on/off switch.

Power Supply: 2x AA Batteries (not included in set). Packing: In Gray or Black pouch with 10 disposable Specula Dia

2.5mm & 4.0 mm
Available All Colours:





06

Otoscope
Ophthalmoscope

METAL OTOSCOPE



TGIND:5200

Fiber Optic
(Metal Otoscope)

Equipped with an easy to use Bayonet Connecting system with high quality Professional metal head and chrome plated knurled battery handle making the surface a secure grip in the hand with extra bright light of 2.5v LED bulb illumination provided through Fiber Optic system.

Three times high quality magnifying lens.

Requires: 2xC Batteries (Not included with set)

Paching: In Black soft zipper pouch with 10 Specula Having Dia 2.5mm & 4.0mm



TGIND:5201

Conventional
(Metal Otoscope)

Equipped with an easy to use Bayonet Connecting system with high quality Professional metal head and chrome plated knurled battery handle making the surface a secure grip in the hand with extra bright light of 2.5v LED bulb illumination provided through Fiber Optic system.

Three times high quality magnifying lens.

Requires: 2xC Batteries (Not included with set)

Paching: In Black soft zipper pouch with 10 Specula Having Dia 2.5mm & 4.0mm



METAL OTOSCOPE



TGIND:5202
Conventional
(Parker Metal Otoscope)

High Quality Professional Parker Otoscope with extra bright light 2.5v LED illumination.
Three times high quality magnifying lens.
Available in Satin and Polish Finish.
Requires: 2xC Batteries (Not included with set)
Paching: In Black soft zipper pouch with 10 Specula Having Dia 2.5mm & 4.0mm



TGIND:5203
Conventional
(Metal Otoscope)

Standard Metal Otoscope with bayonet connecting head fitting.
Ideal for pneumatic tympanic membrane test and every day ENT use.
2.5v LED white Light.
Available in Satin and Polish Finish.
Requires: 2xC Batteries (Not included with set)
Paching: In Black soft zipper pouch with 10 Specula Having Dia 2.5mm & 4.0mm

METAL OTOSCOPE



TGIND:5204

Conventional
(Parker Metal Otoscope)

High Quality Professional Parker Otoscope with screw fitting illumination light. 2.5v LED Metal Handle & Plastic bead CPP. Rheostat dimmer
Requires: 2xC Batteries (Not included with set)



TGIND:5205

Conventional
(Metal Otoscope)

2.5v standard illumination LED bulb. Chrome Plated Handle and head 4x magnifying glasses. Reusable speculum size 0 2.5 3.5 and 4mm include. Connecting for pneumatic Otoscope. Battery handle with rheostat for adjustment of illumination. Supply: German bulb technology for super view with 3 reusable ear specula.
Requires: 2xC Batteries (Not included with set)



METAL OTOSCOPE



TGIND:5206

Conventional
(Veterinary Metal Otoscope)

Open system with rotatable lens and optional speculum for convenient procedures. 2.5v Otoscope.

Requires: 2xC Batteries (Not included with set)

Paching: In Black soft zipper pouch with 3 VET reusable Specula Dia 4.0, 5.0, 7.0mm



TGIND:5207

Conventional
(Veterinary Metal Otoscope)

Code: TI-017-417

Having special focal distance with two fold magnification. Bayonet closure. 2.5v Otoscope.

Requires: 2xC Batteries (Not included with set)

Paching: In Black soft zipper pouch with 3 VET reusable Specula Dia 4.0, 5.0, 7.0mm



DERMATOSCOPE



TGIND:5213

(Dermatoscope)

Bayonet locking system for easy engagement and disengagement from the handle. 2.5v LED illumination of the examination field. High resolution achromatic lenses system with 10x magnification. Adjustable focusing Eyepiece for a large correction range. Lens without scale. Soft eyepiece protection.

Available in Polish and Sand finish.

Power Supply: 2 x C Batteries (not included in set).

Packing: Hard plastic case packing.



TGIND:5214

(Dermatoscope)

Modern pocket Dermatoscope in a compact design with high-quality optics. 10x magnification large, high-resolution image. Choice of two skin-friendly and sterilizable contact plates: with scale from 0-10mm for exact measurement of pigmented skin lesions, without scale. Simple lamp change on the underside of the instrument's head. Enhanced LED Technology for more light. Attachment clip with integrated On/Off switch. High-quality 2 color handle. Power Supply: 2 x AA Batteries (not included in set). Packing: Hard Plastic case packing.

OTOSCOPE LARYNGOSCOPE LED BLUBS

TGIND:5198

LED Bright Bulb For Pediatric & Adult Laryngoscope
50 Gold & 50 Silver With Transparent Box



TGIND:5199

LED Bright Bulb For Otoscope Ophthalmoscope
100Pcs Silver With Transparent Box



DIAGNOSTIC PEN TORCH

Key Features of an Examination Pen Torch

Compact Size: Pen torches are lightweight and easy to carry, making them convenient for use during check-ups or emergencies. **Focused Beam:** The torch emits a concentrated beam of light, which offers precise illumination of areas like the eyes, throat, and other parts of the body.

Adjustable Light: In addition, some models allow users to adjust brightness or focus to meet different examination needs. **Power Source:** Most pen torches use batteries, though rechargeable versions are also available.

Durability: Made from durable materials such as aluminum or plastic, these pen torches can withstand regular use, ensuring their longevity.

Common Uses of an Examination Pen Torch

Eye Exams: Optometrists use pen torches to check pupil reactions, assess retina health, and examine other aspects of the eye.

Throat and Mouth Examinations: These torches help illuminate the back of the throat or the inside of the mouth for better visibility during physical exams.

Neurological Exams: Doctors use pen torches to assess pupil reactions to light, which is a common neurological test.

TGIND:5216

TGIND:5217



www.tgind.com



TGIND Is Adding To It's Tradition Of Providing Quality Diagnostic Instruments



+92-52-3257151
+92-52-3257151
www.tgind.com

TAHIR GROUP OF INDUSTRIES
TAHIR GROUP OF INDUSTRIES
TAHIR GROUP OF INDUSTRIES

info@tgind.com/tgindsale@gmail.com
New Miana Pura East, Roras Road
P.O.S.I.E Sialkot 51310 Pakistan

